



1

Abstract

Recent advancements in non-viral gene editing technologies offer precise genome engineering at significantly reduced costs and biosafety risks compared to viruses. CRISPR is a programmable genome targeting system, which relies on a donor DNA as a template for homology-directed repair (HDR) to introduce exogenous sequences for cell and gene therapy. Therefore, as viruses are increasingly being replaced by safer and cheaper technologies, there is a growing need for reliable and effective vectors for targeted gene insertions. Plasmid DNA (pDNA), which has traditionally been used as the incumbent template for gene-length knock-ins, however, is associated with high toxicities, poor immunogenicity profiles and low HDR efficiencies, all of which has served to limit the success of non-viral gene therapy.

Touchlight's doggybone DNA (dbDNA™) is a fully synthetic, GMP-grade, covalently closed, linear vector, which directly addresses these challenges. The enzymatic manufacture eliminates bacterial sequences and produces a minimal vector of higher purity and lower immunogenicity compared to fermentation-derived DNA. Furthermore, the platform offers unparalleled scalability with reduced timelines overcoming pDNA manufacture bottlenecks.


dbDNA™ has been robustly evaluated as an HDR template in ex vivo CRISPR gene editing platforms internally and by gene therapy industry leaders and is shown to maintain significantly better cell viability compared to pDNA. Improved cell survival leads to accelerated recovery and expansion in culture to produce higher edited primary immune cell counts. With its reduced size, dbDNA™ has a copy-number advantage over pDNA, allowing improved CRISPR knock-in rates at a lower template concentration, further minimising the DNA bioburden. Importantly, the improved toxicity and immunogenicity profiles lead to better and more reliable knock-in efficiencies across primary T cells isolated from different blood donors, reducing undesired donor-to-donor variability. This altogether makes dbDNA™ highly suitable for autologous ex vivo gene and cell therapies, where cell numbers are typically limited and genome editing outcomes inconsistent.

dbDNA™'s improved performance in cells, alongside Touchlight's robust and rapidly scalable manufacturing technology, offers a reliable and effective alternative to pDNA as a knock-in template for CRISPR and, thus, promises to advance non-viral gene therapy.


2

Results


dbDNA's improved viability compared to plasmid allows for increased HDR template concentrations to improve knock-in in primary T cells



dbDNA's improved viability compared to plasmid allows for faster cell recovery and expansion



dbDNA has improved consistency across blood donors



2

Introduction

Doggybone DNA (as a template for gene editing)



Touchlight's doggybone DNA (dbDNA™)

Linear, covalently closed DNA vector

GMP production process in a week

No bacterial sequences

Lower toxicity & immunogenicity




Scan for poster #360 on Touchlight's dbDNA™ technology


- Minimal vector: copy number advantage over plasmid DNA for further reduced toxicity
- Supporting larger cargo capacity > 20kb – more complex and multiple genes
- High purity, low endotoxin, closed ended dsDNA – better editing efficiency
- Avoiding viral vector supply chain – cost, complexity & timeline
- Excellent template for (CRISPR) targeted genome insertions
- A bulk of CRISPR editing data in primary T cells demonstrating dbDNA is superior to other non-viral vectors

dbDNA preserves the physiological T cell profile better than plasmid


CRISPR edited T cell characterisation




Addition of nuclease target sites to dbDNA may improve HDR yields



dbDNA outperforms linear dsDNA in primary T cells



dbDNA & dbDNA+CTS outperform plasmid in primary B cells



4

FAQs & More Touchlight technology

- dbDNA works well for stem cell editing:
 - CRISPR editing data
 - dbDNA has improved viability over pDNA in stem cells
- dbDNA works with a variety of nuclease targeting systems:
 - Cas9 (T cells, B cells, iPSCs)
 - Cas12a (T cells, iPSCs)
 - STAR-CRISPR (T cells)
- dbDNA has shown lower immunogenicity
 - dbDNA formulates well in LNP compared to pDNA:
 - Smaller size
 - Improved DNA cargo
 - Lower level of immunogenicity markers (TLR9 pathway) in iPSC generated with dbDNA compared to those generated using pDNA

dbDNA as a vector for Lentivirus Production – scan for poster #233



mbDNA: Touchlight's circular ssDNA HDR template

